



# Recombinant Streptomyces coelicolor Probable glycerol uptake facilitator protein (glpF)

<b>Product Code</b>	CSB-CF324786FOB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P19255
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Streptomyces coelicolor (strain ATCC BAA-471 / A3(2) / M145)
<b>Sequence</b>	MSSSDIFIGETIGTALLILLGGGVCAAVTLKASKARNAGWLAIAFGWGFAVMTA AYISGPLSGAHLNPAVTVGIAIKDGDWSNTPTYFAGQLLGAMIGAVLVVVAYY GQFQAHLTDREIVGGPGAQDTTAKSVEAQEKGAGPVLGVFSTGPEIRHTVQN LATEIIGTFVLLLAILTQQLNDEGNLGLILGALITGFVVVSIGLSLGGPTGYAINPV RDLGPRIVHALLPLPNKGGSDWSYAWIPVVGPLIGGAIAGGVYNVAFA
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	glpF
<b>Protein Names</b>	Recommended name: Probable glycerol uptake facilitator protein
<b>Expression Region</b>	1-264
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.