



Recombinant Zea mays Aquaporin TIP3-2 (TIP3-2)

Product Code	CSB-CF860166ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9AT75
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MSTATGVRAGRFRFTVGRSEDATHPDTIRAAISEFIATAIFVFAAEGSVLSLGKM YHDHSTISTAGGLVAVALAHALGLAVAVAVAVNVSSGGHVNPVTFGALVGGRRV SLVRVLYWAAQLLGAVAATLLLRLATGGARPPGFALASGVGDGHAVLLEAVM TFGFVYAYYATVVDPKRGHLGTIAPLAVGFLLGANVLAGGPFDFGAGMNPVRF GPALVGWRWRHHWVYWLGPFLGAGLAGLVYEYLLIPPADAVPHTHQPLAPE DY
Source	in vitro E.coli expression system
Target Names	TIP3-2
Protein Names	Recommended name: Aquaporin TIP3-2 Alternative name(s): Tonoplast intrinsic protein 3-2 ZmTIP3-2 ZmTIP3;2
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.