



Recombinant Escherichia coli Leader peptidase pppA (pppA)

Product Code	CSB-CF671607ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q46836
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MLFDVFQQYPTAMPVLATVGGLIIGSFLNVVIWRYPIMLRQQMAEFHGEMSSA QSKISLA LPRSHCPHCQQQTIRIRDNIPLFSWLMLKGRCDQCQAKISKRYPLVELLTALAFLL ASLVW PESGWGLAVMILSAWLIAASVIDLDHQWLDPVFTQGVLWTGLIAAWAQQSPLT LQDAVTG VLVGFITFYSLRWIAGIVLRKEALGMGDVLLFAALGGWVGALSLPNVALIASCC GLIYAV ITRKGSTTLPGPCLSLGGIATLYLQALF
Source	in vitro E.coli expression system
Target Names	pppA
Protein Names	Recommended name: Leader peptidase pppA
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.