



# Recombinant Salmonella typhimurium Nitrite transporter NirC (nirC)

<b>Product Code</b>	CSB-CF326253SXB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25926
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
<b>Sequence</b>	MFTDTINKCAANAARIARLSANNPLGFVSSAMAGAYVGLGIIIFTLGNLLDPS VRPLVMGATFGIALTLVIIAGSELFTGHTMFLTLGVKAGTISHGQMWAILPQTWL GNLVGSVVFVALLYSWGGSLLPVDTSIVHSVALAKTTAPATVLFFKGALCNWL VCLAIWMAIRTEGTAKFLAIWWCLLAFIASGYEHSVANMTLFFALSDFGHHSDA YTLAGIGHNLLWVTLGNTLSGVVFMGLGYWYATPKSERPAPAKINQPEAAANN
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	nirC
<b>Protein Names</b>	Recommended name: Probable nitrite transporter
<b>Expression Region</b>	1-269
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.