



Recombinant Pig Aquaporin-1 (AQP1)

Product Code	CSB-CF744142PI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PQZ1
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MASEFKKKIFWRAVVAEFLAMTLFIFISIGSALGFQYPVRNNQTSGAAQDNVKV SLAFGLSIATLAQSVGHISGAHLNPAVTLGLLLSCQISVLRVAMYIIAQCVGAIWA TAILSGITSSLPGNSLGLNSLAPGVDSGQGLGIEIIGTLQLVLCVLATTD RRRRDL GGSAPLAIGFSVALGHLLAIDYTGCGINPARSFGSAVITHNFDHWVFWVGPFI GGALAVLIYDFILAPRSSDLTDRVKVWTSGQVEEYDLGDGDDINSRVEMKPK
Source	in vitro E.coli expression system
Target Names	AQP1
Protein Names	Recommended name: Aquaporin-1 Short name= AQP-1 Alternative name(s): Aquaporin-CHIP Water channel protein for red blood cells and kidney proximal tubule
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.