



# Recombinant Glycine max Nodulin-26

<b>Product Code</b>	CSB-CF362572GGV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P08995
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Glycine max (Soybean) (Glycine hispida)
<b>Sequence</b>	MADYSAGTESQEVVNVTKNTSETIQRSDSLVSVPFLLQKLVAEAVGTYFLIFAG CASLVNENYYNMITFPGIAIVWGLVLTVLVYTVGHISGGHFNPVAVTIAFASTRR FPLIQVPAYVVAQLLGSILASGTLRLLFMGNHDQFSGTVPNGTNLQAFVFEFIM TFFLMFVICGVATDNRAVGEFAGIAIGSTLLLNVIIIGGPVTGASMNPARSLGPAF VHGEYEGIWIYLLAPVVGAIAGAWVYNIVRYTDKPLSETTKSASFLKGRAASK
<b>Source</b>	in vitro E.coli expression system
<b>Protein Names</b>	Recommended name: Nodulin-26 Short name= N-26
<b>Expression Region</b>	1-271
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.