



Recombinant *Neosartorya fumigata* Protein alcS (alcS)

Product Code	CSB-CF534019NGT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0YBR5
Product Type	Transmembrane Protein
Immunogen Species	<i>Neosartorya fumigata</i> (strain CEA10 / CBS 144.89 / FGSC A1163) (<i>Aspergillus fumigatus</i>)
Sequence	MDTEQGLKNHTAKTSPHDETAMASLTTIPTSVTLTSAEQFEKLYLSPLTQRQGM LSKQMGN PTPLALGGFVITTTPLSCCLMGWRGATGSGIAFTGPIIFLGGLLVLTSILEFILG NTFP CVVFGTIGAFWFAFGCTMTPAFNAAAPFSTSATDTVAGLSSPDFLNTYAFLFIW MGVLML IFLACATRTRNAVYVAIFTTTLVFGFLSGAYWRLAVADALVGNRLVVAAGACLF VASMLG FYLLVAQLFDSVGLPVRLPVGDLRSRFDWRRAR
Source	in vitro <i>E.coli</i> expression system
Target Names	alcS
Protein Names	Recommended name: Protein alcS
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.