



Recombinant Bacillus halodurans Uncharacterized protein BH0177 (BH0177)

Product Code	CSB-CF885144BQS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9KGC7
Product Type	Transmembrane Protein
Immunogen Species	Bacillus halodurans (strain ATCC BAA-125 / DSM 18197 / FERM 7344 / JCM 9153 / C-125)
Sequence	MKKRRSNLELLLQCFYRKKVYLSVFIVLALYCVNLLIENASSEKNVYDLLMQLTE FHPLT YTIIPTYLVLTAHFSLGKMHYLAFRCKDKRQWYNLNVSCIAIVTTGYSVLIAFI MLMQ SLFVFRFENKWSTFAVDYYTYHATFLMNYSPLVYSIATLLLLWLLLFLGLLFYVI FIWT KSPLVSLLFVLLNIMNAAVTLGKIDTLTPVFFTDREVSIQYVYKFDLNQDSFPYSI FVY WIMLIAVIYLIGWLVIQRVDFESEKGEKHHAS
Source	in vitro E.coli expression system
Target Names	BH0177
Protein Names	Recommended name: Uncharacterized protein BH0177
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.