



Recombinant Acidovorax ebreus Undecaprenyl-diphosphatase (uppP)

Product Code	CSB-CF496514AWK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9MF16
Product Type	Transmembrane Protein
Immunogen Species	Acidovorax ebreus (strain TPSY) (Diaphorobacter sp. (strain TPSY))
Sequence	<p>MDLLLLLKAAIMGVVEGLTEFLPISSTGHLILAGSLLGFDDAKAKVFDIAIQTGAIF AVI LVYWQRIRATLVALPTERQARRFALNVLIGFLPAVLLGLLLGAIKAHLFTPVVV ASTFI LGGFVILWAERRQQAAVRIHAVDDMTPLDALKVGLVQCLAMVPGTSRSGATII GGMLLGL SRKAATDYSFFLAIPTLIGAGVYSLYKERALLSAADIPLFAVGLVFSFISAWLCVR WLLR YISSHSFVPFAWYRIAFGLVVLVTAWSGLVTWAE</p>
Source	in vitro E.coli expression system
Target Names	uppP
Protein Names	Recommended name: Undecaprenyl-diphosphatase EC= 3.6.1.27 Alternative name(s): Bacitracin resistance protein Undecaprenyl pyrophosphate phosphatase
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.