



Recombinant Escherichia coli Nickel/cobalt efflux system rcnA (rcnA)

Product Code	CSB-CF631130EGW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1R9W6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Sequence	MTEFTTLLQQGNAWFFIPSAILLGALHGLEPGHSKTMMAAFIIAIKGTIKQAVML GLAAT ISHTAVVWLIAFGGMVISKRFTAQSAEPWLQLISAVIIIISTAFWMFWRTWRGER NWLENM HEHDHEHHHHDHEDHHDHGHHHHHEHGEYQDAHARAHANDIKRRFDGREV TNWQILLFGL TGGLIPCPAAITVLLICIQLKALTLGATLVVSFSLGLALTLVTVGVGAAISVQQVAK RWS GFNTLAKRAPYFSSLLIGLVGVYMGVHGFMGIMR
Source	in vitro E.coli expression system
Target Names	rcnA
Protein Names	Recommended name: Nickel/cobalt efflux system rcnA
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.