



Recombinant Spermidine/putrescine transport system permease protein PotB (potB)

Product Code	CSB-CF360279SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFK5
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MIVTIVGWLVLVFLPNLMIIGTSFLTRDDASFVKMVFTLDNYTRLLDPLYFEVLL HSLN MALIATLACLVLGYPFWFLAKLPHKVRPLLLFLLIVPFWTNSLIRIYGLKIFLSTK GYL NEFLLWLGVIDTPIRIMFTPSAVIIGLVYILLPFMVMPLYSSIEKLDKPLLEAARDL GAS KLQTFIRIIIIPLTMPGIIAGCLLVMLPAMGLFYVSDLMGGAKNLLIGNVIKQVFLNI RDW PFGAATSITLTIVMGLMMLLVYWRASRLLNKKVELE
Source	in vitro E.coli expression system
Target Names	potB
Protein Names	Recommended name: Spermidine/putrescine transport system permease protein PotB
Expression Region	1-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.