



Recombinant Caldicellulosiruptor sp. Putative ABC transporter permease protein ORF2

Product Code	CSB-CF334796CAL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40980
Product Type	Transmembrane Protein
Immunogen Species	Caldicellulosiruptor sp. (strain Rt8B.4)
Sequence	MSFIKETKKWYFIFLAVWTLIADVPFLFMLFTSFKTQSELLSGNTWQIPRQPTIG NFSTV LEGNFFTYLKNSVIAVSISVVLILIISSMAAFAFSRFKFALNNLLYSLIAGMAIPIHV T LIPIYVLTNLIKLYDTVFALIGPYVALSLPMSIFILTEFMREIPLELEEAAKIDGCSM FR LYSDILLPLSAPALITVGIYNGTYLWNEFVFALVLTSSPTRRTLPLGIWDFQARY GSDIP AIMAFLTSLPMLIAYIFGQDKIIKGMMAGAVKG
Source	in vitro E.coli expression system
Protein Names	Recommended name: Putative ABC transporter permease protein ORF2
Expression Region	1-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.