



Recombinant Escherichia fergusonii Rhomboid protease glpG (glpG)

Product Code	CSB-CF485065EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LSC6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Sequence	MLMITSFANPRVAQAFVDYMATQGVILTIIQQHNQSDVWLADESLAERVSELA RFLENPA DPRYLAASWQAGHTGSGLYHRRYPFFAALRERAGPVTWVMMIACVVVFIAMQ ILGDQEVN LWLAWPFDPTLKFEFWRYFTHALMHFSLMHILFNLLWWWYLGGAWEKRLGSG KLIVITLI SALLSGYVQQKFSGPWFGGLSGVVYALMGYVWLRGERDPQSGIYLQRGLIIF ALIWIVAG WFDLFGMSMANGAHIAGLAVGLAMAFVDSL NARKRK
Source	in vitro E.coli expression system
Target Names	glpG
Protein Names	Recommended name: Rhomboid protease glpG EC= 3.4.21.105 Alternative name(s): Intramembrane serine protease
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.