



Recombinant Danio rerio Transmembrane protein 41B (tmem41b)

Product Code	CSB-CF023838DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1A5V7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MAKKRAGNRETESSPLVEQEPRPSKETPVPKGAQSPGGASARMSILLLVVIFA CSACVMY LVFRNFPQLSEDEREKIKIPKDMEDAKALGTVLSKYKDTYYTQVLLAYFATYIFL QTFAI PGSIFLSILSGYLYPFPLALFLVCLCSGLGASFCYMLSylvGRPMVYKYLTERA QKWSQQ VDKHREHLINYIIFLRITPFLPNWFINITSPVINVPLGVFFLGTFLGVAPPSFVAINA GT TLYKLTTAGEAVSWNSLLVLGVLAVVSILPVCFQKKLQKLE
Source	in vitro E.coli expression system
Target Names	tmem41b
Protein Names	Recommended name: Transmembrane protein 41B
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.