



Recombinant FIT family protein ZK265.9 (ZK265.9)

Product Code	CSB-CF688566CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5CZ37
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSTRRSSTRADSTTKRPASPSTPNAALGIFVAIARQILFIDARKVALFYLAFTVLSFI ESRIELDSTYYLVQKHSVLNQYGVKMGWFWTLVIVGPFIFWSSKAHNRRDRD QPIVDVCR LGVGTACWYFSVQFFHKVLALTSMDKGRTRLTRAQCSEKEGVWTPGYDISGH CFLMIYSI LIITEEAIAYRHYQQVTDVAHQMDGDREEHDLTRCIQYFFVAMFLHAFWFKQ IISVL YYHIFIEEILGAVAAVVCWFVTYRMLYPAGFLASPIRRTVGRK
Source	in vitro E.coli expression system
Target Names	fitm-2
Protein Names	Recommended name: FIT family protein ZK265.9
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.