



Recombinant Candida albicans Rhomboid protein 2 (RBD2)

| | |
|--------------------------|--|
| Product Code | CSB-CF685495CZD |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q5AH12 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) |
| Sequence | <pre>MINFDEFPPKYDPSLRPPALTTGLIVFTFMLCVIKSVTSYTFESLILYPRAPLDLN LNSI SLYSLFHVNFFHWICNIFTLATPLAVFETRNGTIHTGVTLNLLTVIAALQYCIVGLI FYP NTGVIGLSGIAFSLMSY MAYHESKFRPIMHTFHLSNSLEIKLYTLYVPFVVAIVF MILFP SSSLPGHLFGITTYGILLSYGYIDKLYPPSKVITTENKLSSLINFLIVTFYKEEES LV TRGSVGYKPLFNQDIEHGAANAASHGSSAFVGETRVLGTRESTV</pre> |
| Source | in vitro E.coli expression system |
| Target Names | RBD2 |
| Protein Names | Recommended name: Rhomboid protein 2 EC= 3.4.21.- |
| Expression Region | 1-284 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |