



Recombinant *Oryza sativa* subsp. japonica ASC1-like protein 3 (Os03g0264000, LOC_Os03g15750)

Product Code	CSB-CF802917OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q84QC0
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MAIRGPEASSFFPLTLVFSVGFFCARFFLDRLVYKPLAAYCFSSKASKLMNDEV RQAKIV KFSESIWKLTYYGSVQAWVLLIIKQEPWSLDTMQYFEGWPNQYMTSSLMLFY MCQCGFYI YSIFALVAWETRKRKDFAVMMSHHVTSILIGYAYLTGFFRIGTIILALHDASDVFL ETAK LCKYTEKELGASLFFGLFALSWLLLRLIYPPFWIIKTSSYQSIISLRKLEKFPPTLY YIF NTMLLTLLVFHIYWWKLIICLMIMKQLNKNKGQVGEDVRSDEDEE
Source	in vitro E.coli expression system
Target Names	Os03g0264000
Protein Names	Recommended name: ASC1-like protein 3 Alternative name(s): Alternaria stem canker resistance-like protein 3
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.