



# Recombinant *Vicia faba* Aquaporin PIP1.1 (PIP1.1)

<b>Product Code</b>	CSB-CF351097VFJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P61838
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Vicia faba</i> (Broad bean) ( <i>Faba vulgaris</i> )
<b>Sequence</b>	MEGKEEDVRVGANKFPERQPIGTSAQSDKDYKEPPPAPFFEPGELSSWSFW RAGIAEFIATFLFLYITVLTVMGVKRSNMCA SVGIQGIAWAFGGMIFALVYCTA GISGGHINPAVTFGLFLARKLSLTRALYYIVMQCLGAICGAGVVKGFQPKQYQA LGGGANTVAHGTYTKGSLGAEIIGTFVLVYTVFSATDAKRNARDSHVPILAPLPI GFAVFLVHLATIPITGTGINPARSLGAIIYNKDHSWDDHWVFWGPFIGAALAA LYHVVVIRAIPFKSRS
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	PIP1.1
<b>Protein Names</b>	Recommended name: Aquaporin PIP1.1 Alternative name(s): Plasma membrane aquaporin 1 Short name= Aquaporin 1 Plasma membrane intrinsic protein 1a Short name= PIP1a
<b>Expression Region</b>	1-286
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.