



Recombinant Bovine Lipid phosphate phosphohydrolase 2 (PPAP2C)

Product Code	CSB-CF642237BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2HJ61
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MERRWVFLLDVLCVLAALPCAILTFVNTPTYKRGFYCGDDSI RYPYRPDTITH GLMAGV IITATVILVSAGEAYLVYTDRLYSRSDFN NYLAALYKVVGTFLFGAAVSQSLTDL AKYMT GRLRPNFLAVCDPDWSRVNCSAYVQVEVCRGSSANVTESRLSFYSGHSSFG MYCMVFLAL YVQARLCWKWARLLRPTVQFFLVAFALYVGYTRVSDHKHHWSDVLVGLLQGA LVA SLTVR YISDFFKARPPQHCPEEEDLERKPSLSLTLALGETDCNHYGYPVSSS
Source	in vitro E.coli expression system
Target Names	PLPP2
Protein Names	Recommended name: Lipid phosphate phosphohydrolase 2 EC= 3.1.3.4 Alternative name(s): PAP2-gamma Short name= PAP2-G Phosphatidate phosphohydrolase type 2c Phosphatidic acid phosphatase 2c Short name= PAP-2c Short n
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.