



Recombinant Human 3-hydroxyacyl-CoA dehydratase 1 (PTPLA)

Product Code	CSB-CF019018HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0YJ81
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGRLTEAAAAGSGSRAAGWAGSPPTLLPLSPTSPRCAATMASSDEDGTNGG ASEAGEDRE APGERRRRLGVLATAWLTFYDIAMTAGWLVLAIAMVRFYMEKGTHRGLYKSIQK TLKFFQT FALLEIVHCLIGIVPTSVIVTGVQVSSRIFMVWLITHSIKPIQNEESVVLFLVAWTV TEI TRYSFYTFSLLDHLPYFIKWARYNFFIILYPVGVAGELLTIYAALPHVKKTGMFSI RLPN KYNVSFDYYYFLLITMASYIPLFPQLYFHMLRQRRKVLHGEVIVEKDD
Source	in vitro E.coli expression system
Target Names	HACD1
Protein Names	Recommended name: 3-hydroxyacyl-CoA dehydratase 1 Short name= HACD1 EC= 4.2.1.- Alternative name(s): Cementum attachment protein Protein-tyrosine phosphatase-like member A
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.