



Recombinant Escherichia coli O127:H6 UPF0761 membrane protein yihY (yihY)

Product Code	CSB-CF483064EOB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7UNK3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Sequence	MLKTIQDKAKHRTRPLWAWLKLLWQRIDEDNMTTLAGNLAYVSLLSLVPLVAV VFALFAA FPMFSDVSIQLRHFIFANFLPATGDVIQRYIEQFVANSNKMTAVGACGLIVTALL LMYSI DSALNTIWRSKRARPKIYSFAVYWMILTLGPLLAGASLAISSYLLSLRWASDLNT VIDNV LRIFPLLLSWISFWLLYSIVPTIRVPNRDAIVGAFVAALLFEAGKKGFFALYITMFPS YQL IYGVLAVIPILFVWVYWTWCIVLLGAEITVTLGEYRKLKQAAEQEEDDEP
Source	in vitro E.coli expression system
Target Names	yihY
Protein Names	Recommended name: UPF0761 membrane protein yihY
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.