



Recombinant Human Kidney mitochondrial carrier protein 1 (SLC25A30)

Product Code	CSB-CF732909HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5SVS4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSALNWKPFVYGGGLASITAECGTFPIDLTKTRLQIQGQTNDAKFKEIRYRGMHLH ALVRIG REEGLKALYSGIAPAMLRQASYGTIKIGTYQSLKRLFIERPEDETLPINVICGILS GVIS STIANPTDVLKIRMQAQSNIIQGGMIGNFMNIYQQEGTRGLWKGVSLSLAQRAAI VVGVEL PVYDITKKHLILSGLMGDTVYTHFLSSFTCGLAGALASNPVDVVRTRMMNQRV LRDGRCS GYTGTLDCLLQTWKNEGFFALYKGFWPWNLRLGPWNIIFVITYEQLKKLDL
Source	in vitro E.coli expression system
Target Names	SLC25A30
Protein Names	Recommended name: Kidney mitochondrial carrier protein 1 Alternative name(s): Solute carrier family 25 member 30
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.