



Recombinant Danio rerio Solute carrier family 25 member 47-A (slc25a47a)

Product Code	CSB-CF626046DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1ECW7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MHFADFLAGSFGGACGVAVGYPLDVKVRIQTQKQFTGIWQCIVLTIRKEGVH GFFKGMF LPITTISMTSSVVFVGTYRNCLQALSYIRKAENTKLDVFMSSGLAGGVAQVSVMS P GDIVKV RLQCQTESRHSVNPYKYSVKPKYSGPIHCLLSICREQGLSGLYRGALPLALRDG PSFATYF LTYHTLCARLTPDGQKEPEWTVVLLSGGVAGMSGWAVGTPMDVIKARLQMD GVRGQRRYR GLLHCLTVTTRTEGLGVFFRSLGINCLRAFPVNMVVFAVYEVSVRVLRSA PLVS
Source	in vitro E.coli expression system
Target Names	slc25a47a
Protein Names	Recommended name: Solute carrier family 25 member 47-A Alternative name(s): Hepatocellular carcinoma down-regulated mitochondrial carrier homolog A
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.