



Recombinant Dictyostelium discoideum Aquaporin-B (aqpB)

Product Code	CSB-CF706896DKK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54WT8
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	MSLKRSDDYQDLEEGIAMEDGGNIKDEEEKPLDPIEEQNKKRWVLIRAVLGELL CTFLFVYVLCATSANFIRLGSPPNPVVGGLSTGFAAVALIYSFADVSGAHFNPA VTFATCVTRKTSITKGLMYVGAQLVGSVLASLILLATFPGNFPGDKNAASAVAIA PSTDANIGNAFLTELVLTFILVYVIFAVAFDVTDNSVKTKVVGKSSSNLTIYTTS GQTKAGFAPIAIGFTLGFLCFLGGSVSGGAFNPARVFGTALVGNNWTRHWMY WIADFLGAGLAGFAQKFFSSTHK
Source	in vitro E.coli expression system
Target Names	aqpB
Protein Names	Recommended name: Aquaporin-B
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.