



Recombinant *Vibrio vulnificus* Maltose transport system permease protein malG (malG)

Product Code	CSB-CF762652VCQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7MFC1
Product Type	Transmembrane Protein
Immunogen Species	<i>Vibrio vulnificus</i> (strain YJ016)
Sequence	MAMVQGKSLKYRVWATHIALWAFLSMIIFPLLMIVAISFREGNFATGSLIPDNPS LEHWK LALGFSVTNADGSVTPPPFPVLTWLWNSVKVAGITSILIVALSTTSAYAFARLRF KGET ILKAMMIFQMFPVAVLALVALYALFDKLGQYIPFLGLNTHGGLIFSYLGGIALHW TIKGY FETIDNSLEEAAALDGATPWQAFRLVLLPLSVPILAVVFILSFIVVGEVPVASLL LSDV NSYTLAVGMQQYLYPQNYLWGDFAAAAVLSALPITIVFLLAQRWLVGGLTAGG VKG
Source	in vitro E.coli expression system
Target Names	malG
Protein Names	Recommended name: Maltose transport system permease protein malG
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.