



Recombinant Arabidopsis thaliana LAG1 longevity assurance homolog 2 (LAG2)

Product Code	CSB-CF873282DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9LJK3
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MESVSSRGGDPVVKPSMEVWHFQIAVYFAFGFFFLRLVLDRYVFQRIALWLLS TGSAPIK LNDAAATRAKIVKCKESLWKLLYYAACDFFVLQVIYHEPWARDIKLYFHGWPNQ ELKLSIK LYYMCQCGFYVYGVAALLAWETRRKDFAVMMSHHVITIILLSYSYLTSSFFRIGAI ILALH DASDVFMETAKIFKYSEKEFGASVCFALFAVSWLLLRLIYFPFWIIRATSIELLDY LDMT SAEGTLMYYSFNTMLLMLLVFHIYWWYLICAMIVRLLKNRGKVGEDIRSDSEDD DD
Source	in vitro E.coli expression system
Target Names	LAG2
Protein Names	Recommended name: LAG1 longevity assurance homolog 2 Short name= LAG1 homolog 2
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.