



Recombinant Peptide transport system permease protein sapC (sapC)

Product Code	CSB-CF313438SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGH7
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MPYDSVYSEKRPPGTLRTAWRKFYSDASAMVGLYGCAGLAVLCIFGGWFAPY GIDQQFLG YQLLPPSWSRYGEVSFFLGTDDLGRDVLRSLLSGAAPT VGGAFVVTLAATICG LVLGTFA GATHGLRSAVLNHILDTLAIPSLLLAIIVVAFAGPSLSHAMFAVWLALLPRMVR SIYSM VHDELEKEYVIAARLDGASTLNILWFAVMPNITAGLVTEITRALSMAILDIAALGF LDLG AQLPSPEWGAMLGDALELIYVAPWTVMLPGAAIMISVLLVNULLGDGVRRAIAG VE
Source	in vitro E.coli expression system
Target Names	sapC
Protein Names	Recommended name: Peptide transport system permease protein sapC
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.