



Recombinant *Saccharomyces cerevisiae* Carrier protein YMC2, mitochondrial (YMC2)

Product Code	CSB-CF336398SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38087
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	VLKDIFAGTIGGIAQVLVGQPFDTTKVRLQTATTRTTTTLEVLRLNVKNEGVFAFY KGALT PLLGVGICVSVQFGVNEAMKRFFQNYNASKNPNMSSQDVDLSRSNTLPLSQY YVCGLTGG VVNSFLASPIEQIRIRLQTQTSNGGDREFKGPWDCIKKLKAQGGLMRGLFPTMI RAGHGL GTYFLVYEALVAREIGTGLTRNEIPPWKLCLFGAFSGTMLWLTVYPLDVVKSIIQ NDDLRL KPKYKNSISYVAKTIYAKEGIRAFFKGFGPTMVRSA PVNGATFLTFELVMRFLG EE
Source	in vitro E.coli expression system
Target Names	YMC2
Protein Names	Recommended name: Carrier protein YMC2, mitochondrial
Expression Region	34-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.