



Recombinant Pongo abelii ADP/ATP translocase 2 (SLC25A5)

Product Code	CSB-CF735671PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R5A1
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	TDAAVSFAKDFLAGGVAAAISKTAVAPIERVKLLLQVQHASKQITADKQYKGIID CVVRI PKEQGVLSFWRGNLANVIRHFPTQALNFAFKDKYKQIFLGGVDKRTQFWRYF AGNLASGG AAGATSLCFVYPLDFARTRLAADVKGAGAEREFRLGDCLVKIYKSDGIKGLY QGFNVSV QGIIYRAAYFGIYDTAKGMLPDPKNTHIVISWMIAQTVTAVAGLTSYPFDTVRR RMMMQ SGRKGTDIMYTGTLDCWRKIARDEGGKAFFKGAWSNVLRGMGGAFALVLYDE IKKYT
Source	in vitro E.coli expression system
Target Names	SLC25A5
Protein Names	Recommended name: ADP/ATP translocase 2 Alternative name(s): ADP,ATP carrier protein 2 Adenine nucleotide translocator 2 Short name= ANT 2 Solute carrier family 25 member 5
Expression Region	2-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.