



# Recombinant Pig ADP/ATP translocase 3 (SLC25A6)

<b>Product Code</b>	CSB-CF747526PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q6QRN9
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Sequence</b>	TEQAISFAKDFLAGGIAAAISKTAVAPIERVKLLLQVQHASKQIAADKQYKGVDC IVRI PKEQGVLSFWRGNLANVIRYFPTQALNFAFKDKYKQIFLGGVDKHTQFWRYFA GNLASGG AAGATSLCFVYPLDFARTRLAADV GKSATEREFKGLGDCLVKITKSDGIRGLYQ GFNVSV QGIIYRAAYFGVYGTA KGMLPDRNTHIVVSWMIAQTVTAVAGVFSYPFDTVR RRMMM Q SGRKGADIMYKGTLD CWRKIFKDEGGKAFFKGAWSNVLRGMGGAFVLVLYDE LKKVI
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SLC25A6
<b>Protein Names</b>	Recommended name: ADP/ATP translocase 3 Alternative name(s): ADP,ATP carrier protein 3 Adenine nucleotide translocator 3 Short name= ANT 3 Solute carrier family 25 member 6
<b>Expression Region</b>	2-298
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.