



Recombinant Bovine ADP/ATP translocase 3 (SLC25A6)

Product Code	CSB-CF021520BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32007
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	TEQAISFAKDFLAGGIAAAISKTAVAPIERVKLLLQVQHASKQIAADKQYKGVDC IVRI PKEQGVLSFWRGNLANVIRYFPTQALNFAFKDKYKQIFLGGVDKRTQFWRYFA GNLASGG AAGATSLCFVYPLDFARTRLAADVKGSGSEREFRGLGDCLVKITKSDGIRGLY QGFNVSV QGIIYRAAYFGIYDTAKGMLPDPKNTHIVVSWMIAQTVTAVAGVVSYPFDTVRR RMMMQ SGRKGADIMYKGTVDCWRKILKDEGGKAFFKGAWSNVLRGMGGAFVLVLYD ELKKVI
Source	in vitro E.coli expression system
Target Names	SLC25A6
Protein Names	Recommended name: ADP/ATP translocase 3 Alternative name(s): ADP,ATP carrier protein 3 ADP,ATP carrier protein, isoform T2 Short name= ANT 2 Adenine nucleotide translocator 3 Short name= ANT 3 Solute carrier family 25 mem
Expression Region	2-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.