



# Recombinant Tachyglossus aculeatus aculeatus ADP/ATP translocase 2 (SLC25A5)

<b>Product Code</b>	CSB-CF021518TKU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q000K2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Tachyglossus aculeatus aculeatus (Southeast Australian short-beaked echidna)
<b>Sequence</b>	TDAAVSFAKDFLAGGVAAAISKTAVAPIERVKLLLQVQHASKQITADKQYKGIM DCVVRI PKEQGVLSFWRGNLANVIRYFPTQALNFAFKDKYKQIFLGGVDKRTQFWRYFA GNLASGG AAGATSLCFVYPLDFARTRLAADVKGAGDAREFKGLGDCLVKITKSDGIRGLY QGFNVSV QGIIYRAAYFGIYDTAKGMLPDPKNTHIFISWMIAQSVTAVAGLTSYPFDTVRR RMMMQ SGRKGSDIMYTGIDCWKKIARDEGSKAFFKGAWSNVLRGMGGAFVLVLYDEI KKYT
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SLC25A5
<b>Protein Names</b>	Recommended name: ADP/ATP translocase 2 Alternative name(s): ADP,ATP carrier protein 2 Adenine nucleotide translocator 2 Short name= ANT 2 Solute carrier family 25 member 5
<b>Expression Region</b>	2-298
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.