



Recombinant Mouse Mitochondrial 2-oxodicarboxylate carrier (Slc25a21)

Product Code	CSB-CF803965MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BZ09
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSASNVSLHETSRQVAAGGSAGLVEICLMHPLDVVKTRFQVQRSVTDPQSY RTVRGSFQ MIFRTEGLFGFYKGIIPPILAETPKRAVKFSTFELYKKFLGYMSLSPGLTFLIAGL GSGL TEAVVVNPFVVKVGLQVNRNLFKEQPSTFAYARQIIKKEGLGFQGLNKGLTAT LGRHGI FNMVYFGFYHNVKNIIPSSKDPTLEFLRKFGIGFVSGTMGSVFNIPFDVAKSRIQ GPQPV PGEIKYRSCFKTMEMIYREEGILALYKGLVPMRLGPGGGVMLLVYEYTYAW LQENW
Source	in vitro E.coli expression system
Target Names	Slc25a21
Protein Names	Recommended name: Mitochondrial 2-oxodicarboxylate carrier Alternative name(s): Solute carrier family 25 member 21
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.