



Recombinant Dipeptide transport system permease protein dppC (dppC)

Product Code	CSB-CF360108EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AEG2
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MSQVTENKVISAPVPM TPLQEFWHYFKRNKGAVVGLVYVVIVLFIAIFANWIAP YNPAEQ FRDALLAPPAWQEGGSMALLGTDDVGRDVL SRLMYGARLSLLVGCLVVVLS LIMGVILG LIAGYFGGLVDNIIMRVVDIMLALPSLLLALVLAIFGPSIGNAALALTFVALPHYV RLT RAAVLVEVNRDYVTASRVAGAGAMRQMFINIFPNCLAPLIVQASLGFSNAILDM AALGFL GMGAQPPTPEWGTMLSDVLQFAQSAWWVTFPGLAILLTVLAFNLMGDGLR DALDPKCLKQ
Source	in vitro E.coli expression system
Target Names	dppC
Protein Names	Recommended name: Dipeptide transport system permease protein dppC
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.