



Recombinant Glutathione transport system permease protein gsiD (gsiD)

Product Code	CSB-CF820219EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8X6V6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MRLFNWRRQAVLNAMPLVKPDQVRTPWHEFWRRFRRQHMAMTAALFVILLIV VAIFARWI APYDAENYFDYDNLNNGPSLQHWFGVDSLGRDIFSRVLVGAQISLAAGVFAVF IGAAIGT LLGLLAGYYEGWWDRLIMRICDVLFAPPGILLAIAVVAVLGSGIANVIIAIVAIFSIP AFA RLVRGNTLVLKQQTFIESARSIGASDMTILLRHILPGTVSSIVVFFTMRIGHTSIISA ASL SFLGLGAQPPTPEWGAMLNEARADMVIAPHVAVFPALAIFLTVLAFNLLGDGLR DALDPK IKG
Source	in vitro E.coli expression system
Target Names	gsiD
Protein Names	Recommended name: Glutathione transport system permease protein gsiD
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.