



Recombinant Escherichia coli Glutathione transport system permease protein gsiD (gsiD)

Product Code	CSB-CF300512ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75799
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MRLFNWRRQAVLNAMPLVKPDQV RTPWHEFWRRFRRQHMAMTAALFVILLIV VAIFARWI APYDAENYFDYDNLNNGPSLQHWFGVDSLGRDIFSRVLVGAQISLAAGVFAVF IGAAIGT LLGLLAGYYEGWWDRLIMRICDVLFAPPGILLAI VVAVLGSGIANVIIA VAIFSIP AFA RLVRGNTLV LKQQT FIESARSIGASDMTVLLRHILPGTVSSIVVFFTM RIGTSIISA ASL SFLGLGAQPPTPEWGAMLNEARADMVIAPHVAVFPALAI FLTVLAFNLLGDGLR DALDPK IKG
Source	in vitro E.coli expression system
Target Names	gsiD
Protein Names	Recommended name: Glutathione transport system permease protein gsiD
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.