



Recombinant Bacillus pseudofirmus Dipeptide transport system permease protein dppC (dppC)

Product Code	CSB-CF309034BRB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P94312
Product Type	Transmembrane Protein
Immunogen Species	Bacillus pseudofirmus (strain OF4)
Sequence	MKTEHAKPLMTPEPNSPPPEDQYTAVSPWREAFKQLKKNKLAIVGLIITSFILIA IFAP LLTSSSYAETNPSNRLQGPSAEHWFGTDDFGRDIFTRIVYGARLSLQVGGFAV TGALIFG TTLGLIAGYYGRWIDMLISRIFDIMLAFPSILLAIIVAILGPSLQNALIAIAIVNVPIF GRLVRSKVISLREEEFIMAACAQGMKNGRIIFHHILPNSLAPIIVQATLGFGTALIE AAA LGFLGLGAQAPMPEWGWKMLSDSRQFIQSAPWTVLFPGFSSIMLVVLGFNMIGD GLRDALDP KMKS
Source	in vitro E.coli expression system
Target Names	dppC
Protein Names	Recommended name: Dipeptide transport system permease protein dppC
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.