



Recombinant Oligopeptide transport system permease protein oppB (oppB)

Product Code	CSB-CF360264SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AFH5
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MLKFILRRCLEAIPTLFILITISFFMMRLAPGSPFTGERTLPPEVMANIEAKYHLN DPIM TQYFSYLKQLAHGDFGPSFKYKDYSVNDLVASSFPVSAKLGAAFFLAVILGVS AGVIAA LKQNTKWDYTVMGLAMTGVVIPSFVVAPLLVMIFAILHWLPGGGWNGGALKF MILPMVA LSLAYIASIARITRGSMIEVLHSNFIRTARAKGLPMRRIILRHALKPALLPVLSYMG PAF VGIITGSMVIETIYGLPGIGQLFVNGALNRDYSLVLSLTILVGALTILFNAIVDVLYA VI DPKIRY
Source	in vitro E.coli expression system
Target Names	oppB
Protein Names	Recommended name: Oligopeptide transport system permease protein oppB
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.