



Recombinant Salmonella typhimurium Glutathione transport system permease protein gsiC (gsiC)

Product Code	CSB-CF838493SXB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8ZQM2
Product Type	Transmembrane Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Sequence	MLNYVLKRLGLIPTLLIVAVLVFLFVHLLPGDPARLIAGPEADAQVIALVRQQLG LDQP LHVQFWHYITHVLQGDFGTSMVSRPVSEEIASRFLPTLWLTITSMIWAFLFGM AIGIAA AVWRNRWPDRLGMTLAVTGISFPAFALGMLLMQIFSVDLGWLPTVGADSWQ HYILPSLTL GAAVASVMARFTRSSFVDVLSYDMRTARAKGVSETWVVLKHGLRNAMIPVV TMMGLQFG FLLGGSIVVEKVFNWPGGLGRLLVDSVDMRDYPVIAEVLLFSLEFILINLVVDVL YAAIN PAIRYK
Source	in vitro E.coli expression system
Target Names	gsiC
Protein Names	Recommended name: Glutathione transport system permease protein gsiC
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.