



Recombinant Rabbit Mitochondrial brown fat uncoupling protein 1 (UCP1)

Product Code	CSB-CF025554RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14271
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MVGTTTTDVPPTMGVKIFSAGVAACLADVITFPLDTAKVRQQIQGEFPITSGIRY KGLG TITTLAKTEGPKLYSGLPAGLQRQISFASLRIGLYDTVQEFFTSGEETPSLGSKI SAGL TTGGVAVFIGQPTEVVKVRLQAQSHLHGLKPRYTGTYNAYRIIATTESLTSLWK GTTPNL LRNVIINCTELVTYDLMKGALVRNEILADDVPCHFVSALIAGFCTLLSSPVDVV KTRFI NSPPGQYASVPNCAMTMFTKEGPTAFFKGFVPSFLRLGSWNVIMFVCFEKLK GELMRSRQ TVDCAT
Source	in vitro E.coli expression system
Target Names	UCP1
Protein Names	Recommended name: Mitochondrial brown fat uncoupling protein 1 Short name= UCP 1 Alternative name(s): Solute carrier family 25 member 7 Thermogenin
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.