



Recombinant Zea mays ADP,ATP carrier protein 1, mitochondrial (ANT1)

Product Code	CSB-CF356337ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04709
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	APAEEKGGKNFMIDFMMGGVSAAVSKTAAAPIERVKLLIQNQDEMIKSGRLSEP YKGIVDC FKRTIKDEGFSSLWRGNTANVIRYFPTQALNFAFKDYFKRLFNFKKDRDGYWK WFAGNLA SGGAAGASSLFFVYSLDYARTRLANDAKAAKGGGERQFNGLVDVYRKTLLKSD GIAGLYRG FNISCVGIIVYRGLYFGLYDSIKPVVLTGNLQDNFFASFALGWLITNGAGLASYP I D T V R RRMMMTSGEAVKYKSSLDAFQQILKKEGPKSLFKGAGANILRAIAGAGVLSGY DQLQILF FGKKYGSGGA
Source	in vitro E.coli expression system
Target Names	ANT1
Protein Names	Recommended name: ADP,ATP carrier protein 1, mitochondrial Alternative name(s): ADP/ATP translocase 1 Adenine nucleotide translocator 1 Short name= ANT 1
Expression Region	78-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.