



Recombinant Arabidopsis thaliana ADP,ATP carrier protein 3, mitochondrial (AAC3)

Product Code	CSB-CF529869DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O49447
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	APSEKTGTGFLIDFLMGGVSAAVSKTAAAPIERVKLLIQNQDEMIKAGRLSEPY KGISDC FARTVKDEGMLALWRGNTANVIRYFPTQALNFAFKDYFKRLFNFKKEKDGW KWFAGNLA SGGAAGASSLLFVYSLDYARTRLANDAKAAKGGQRQFNGMVDVYKKTIASD GIVGLYRG FNISCVGIVVYRGLYFGLYDSLKPVVLVDGLQDSFLASFLLGWGITIGAGLASYP IDTVR RRMMMTSGEAVKYKSSLQAFSQIVKNEGAKSLFKGAGANILRAVAGAGVLAG YDKLQLIV LGKKYGSGGG
Source	in vitro E.coli expression system
Target Names	AAC3
Protein Names	Recommended name: ADP,ATP carrier protein 3, mitochondrial Alternative name(s): ADP/ATP translocase 3 Adenine nucleotide translocator 3 Short name= ANT 3
Expression Region	70-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.