



Recombinant Rat Lipid phosphate phosphohydrolase 3 (Ppap2b)

Product Code	CSB-CF018415RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97544
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MQSYKYDKAIVPESKNGGSPALNNNPRKGGSKRVLLICLDLFLFMAALPFLIIE TSTIK PYRRGFYCNDESIKYPLKVSETINDAVLCAVGIVIAILRIITGEFYRIYYLKEKSRS TIQ NPYVAALYKQVGCFLFGCAISQSFTDIAKVSIGRLRPHFLSVCDPDFSQINCSE GYIQNY RCRGEDSKVQEARKSFFSGHASFSMFTMLYLVLVYLQARFTWRGARLLRPLLQ FTLLMMAF YTGLSRVSDYKHHPDVLGFAQGALVACCIVFFVSDLFKTKTTLSPAPAIRR EILSPV DIMDRSNHHNMV</p>
Source	in vitro E.coli expression system
Target Names	Plpp3
Protein Names	<p>Recommended name: Lipid phosphate phosphohydrolase 3 EC= 3.1.3.4 Alternative name(s): Differentially expressed in rat intestine 42 Short name= Dri42 PAP2-beta Phosphatidate phosphohydrolase type 2b Phosphatidic acid phosp</p>
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>