



Recombinant Zinc transporter zitB (zitB)

Product Code	CSB-CF769520SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q83SA2
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MAHSHSHTSSHLPEDNNARRLLYAFGVTAGFMLVEVVGGFLSGSLALLADAG HMLTDTAA LLFALLAVQFSRRPPTIRHTFGWLRLTTLAAFVNAIALVVITILIVWEAIERFRTPR PVE GGMMMAIAVAGLLANILSFWLLHHGSEEKNLNVRAAALHVLGDLLGSGVGAIIAA LIIIWT GWTPADPILSILVSLVLRSAWRLLKDSVNELLEGAPVSLDIAELKRRMCREIPE VRNVH HVHVWMVGEKPVMTLHVQVIPPHDHDALLDQIQHYLMDHYQIEHTTIQMEYQP CHRPDCH LNEGVS GHSHHHH
Source	in vitro E.coli expression system
Target Names	zitB
Protein Names	Recommended name: Zinc transporter zitB
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.