



Recombinant *Macaca fascicularis* ADP/ATP translocase 4 (SLC25A31)

Product Code	CSB-CF700184MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4R8M0
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MHREPPKKKAEKRLFDASSFGKDLLAGGVAAAVSKTAVAPIERVKLLLQVQAS SKQISPE ARYKGMVDCLVRIPREQGFFSFWRGNLANVIRYFPTQALNFAFKDKYKQLFMS GVNKEKQ FWRWFLANLASGGAAGATSLCVVYPLDFARTRLGVDIGKGPEERQFKGLGDC IMKIAKSD GIAGLYQGFVSVQGIIVYRASYFGAYDTVKGLLPKPKKTPFLVSFFIAQVVTTCC SGILS YPFDTVRRRMMMQSGEAKRQYKGTLDLDFVKIYQHEGINSFFRGAFSNVLVGRGT GGALVVLVLDYKIKEFFHIDIGGR
Source	in vitro E.coli expression system
Target Names	SLC25A31
Protein Names	Recommended name: ADP/ATP translocase 4 Alternative name(s): ADP,ATP carrier protein 4 Adenine nucleotide translocator 4 Short name= ANT 4 Solute carrier family 25 member 31
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.