



Recombinant Phosphate transport system permease protein pstC (pstC)

Product Code	CSB-CF314724SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGI0
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MAATKPAFNPPGKKGDIIFSVLVKLAALIVLLMLGGIIVSLIISSWPSIQKFGLAFL WTK EWDAPNDIYGALVPIYGTLVTSFIALLIAPVVSFGIALFLTELAPGWLKRPLGIAIE LLA AIPSIVYGMWGLFIFAPLFAVYFQEPVGNIMSNIPVIGALFSGPAFGIGILAAGVIL AIM IIPYIAAVMRDVF EQTPVMMKESAYGIGCTTWEVIWRIVLPFTKNGVIGGIMLGL GRALG ETMAVTFIIGNTYQLDSASLYMPGNSITSALANEF AEAESGLHVAALMELGLILF VITFI VLAASKFMIMRLAKNEGAR
Source	in vitro E.coli expression system
Target Names	pstC
Protein Names	Recommended name: Phosphate transport system permease protein pstC
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.