



# Recombinant Mouse ADP/ATP translocase 4 (Slc25a31)

<b>Product Code</b>	CSB-CF668906MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q3V132
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	MSNESSKKQSSKKALFDPVSFSKDLLAGGVAAAVSKTAVAPIERVKLLLQVQA SSKQISP EARYKGMLDCLVRIPREQGFLSYWRGNLANVIRYFPTQALNFAFKDKYKELFM SGVNKEK QFWRWFLANLASGGAAGATSLCVVYPLDFARTRLGVDIGKGPEQRQFTGLGD CIMKIAKS DGLIGLYQGFGVSVQGIIVYRASYFGAYDTVKGLLPKPKETPFLVSFIIAQIVTTC SGIL SYPFDTVRRRMMMMSGESDRQYKGTIDCFLKIYRHEGVPAFFRGAFSNILRGT GGALVLV LYDKIKEFLNIDVGGSSSGD
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Slc25a31
<b>Protein Names</b>	Recommended name: ADP/ATP translocase 4 Alternative name(s): ADP,ATP carrier protein 4 Adenine nucleotide translocator 4 Short name= ANT 4 Solute carrier family 25 member 31 Sperm flagellar energy carrier protein
<b>Expression Region</b>	1-320
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.