



Recombinant Ribose transport system permease protein rbsC (rbsC)

Product Code	CSB-CF314191EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AGI3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MTTQTVSGRRYFTKAWLMEQKSLIALLVLIIVSTLSPNFFTINNLFNILQQTSVN AIMA VGMTLVILTSGIDLSVGSLLALTGAVAASIVGIEVNALVAVAAALALGAAIGAVTG VIVA KGRVQAFIATLVMMLLLRGVTMVYTNGSPVNTGFTENADLFGWFGIGRPLGVP TPVWIMG IVFLAAWYMLHHTRLGRIYALGGNEAATRLSGINVNKIKIIVYSLCGLLASLAGII EVA RLSSAQPTAGTGYELDAIAAVVLGGTSLAGGKGRIVGTLIGALILGFLNNGLNLL GVSSY YQMIVKAVVILLAVLVDNKKQ
Source	in vitro E.coli expression system
Target Names	rbsC
Protein Names	Recommended name: Ribose transport system permease protein rbsC
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.