



Recombinant Protein translocase subunit SecF (secF)

Product Code	CSB-CF313430EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AG95
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MAQEYTVSQLNHGRKVYDFMRWDYWAFGISGLLLIAAIVIMGVRGFNWGLDF TGGTVIEI TLEKPAEIDVMRDALQKAGFEEMQLQNFSSHDIMVRMPPAEGETGGQVLGS QVLKVINE STNQNAAVKRIEFVGPSVGADLAQTGAMALMAALLSILVYVGFRFEWRLAAGV VIALAHD VIITLGILSLFHIEIDLTIVASLMSVIGYSLNDSIVVSDRIENFRKIRRGTPYEIFNV S LTQTLHRTLITSGTTLMLMLYLFGGPVLEGFSLTMLIGVSIGTASSIYVASALAL KLG MKREHMLQQKVEKEGADQPSILP
Source	in vitro E.coli expression system
Target Names	secF
Protein Names	Recommended name: Protein translocase subunit SecF
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.