



Recombinant Schizosaccharomyces pombe Uncharacterized mitochondrial carrier C17H9.08 (SPAC17H9.08)

Product Code	CSB-CF515587SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13805
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MPNSTAQYPEKDSWEFLVKSGIAGGTAGCVAKSVVAPLDRVKILYQTNHASYR GYAYSRH GLYKAIKHIYHVYGLHGLYQGHTATLYRVFPYAGIKFVAYEQVRRVLIRDPEHE THARRF LSGSLAGTCSVFFTYPLELIRVRLAYITNTGKNPTLTQVTKDIFHERDFLCNKKY PGLSR LSKICNFYRGFSVTLTGIFPYAGMSFLAYDLATDFFHKQKIDEWVSTKKSDKKL KTWPEL LCGAFAGVCGQTVSYPFEVCRRKMQIGGIRKNKSFLRLKQVVQTTYKEAGMR GFFVGLTI GYIKVIPMVSTSFVYNHASKALLGID
Source	in vitro E.coli expression system
Target Names	SPAC17H9.08
Protein Names	Recommended name: Uncharacterized mitochondrial carrier C17H9.08
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.